

DERMAMAXX™



The Complete Dermamaxx Bundle is self-contained and required no computer or other devices. Just plug it in and use it. No programming, no PC, no hassle. The high resolution color flat panel screen display and portable hand held camera creates amazing opportunities to see problem areas close up.

With fully adjustable magnification, the self-illuminated camera is cool to the touch, easy to use, and is one of the finest dermatology imaging solutions available.

The Dermamaxx is a complete solution for diagnosing animal dermal conditions. With a micro-camera capable of 100x magnification, the Dermamaxx enables you to get in close to the problem and give a proper diagnosis.

www.dermamaxx.com

1-949-428-3976
1(800)400-7972 (U.S. & Canada Only)

8697 Research Drive
Irvine, California, USA 92618

Micro-Camera

- Lightweight
- Vibrant clarity
- Portable
- Self-contained lighting
- Full-bright images
- Easy-to-use

Integrated LCD

- 17.1" TFT LCD
- 1280 x 1024 Resolution
- 16.7 Million Colors
- Full Brightness
- High Contrast

Special Options

- Mobile Kart
- Image Capture & Store
- Wall Mount
- Poll Mount



 **Ampronix**
IMAGING TECHNOLOGY