

CHLORIDE

PCm Series Medical-Grade Power Conditioner

Table 4. Specifications for 120 VA - 1000 VA Models

Model Number*	Load (VA)	Nominal Input (Volts)	Output (Volts)	Max. Load (Amps)	Frequency (Hz)	Connectors: Input / Output**	Dimensions (HxWxD) in (cm)	Ship Weight lbs (kg)
PCM120A	120	120	120	1.0	60	5-15P HG / (2) 5-15R HG	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	9.5 (4.3)
PCM180A	180	120	120	1.5	60	5-15P HG / (2) 5-15R HG	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	9.5 (4.3)
PCM180I	180	200 - 240	200 - 240	0.8	50/60	IEC320 / (4) IEC320	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	12 (5.4)
PCM240A	240	120	120	2.0	60	5-15P HG / (4) 5-15R HG	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	11.5 (5.2)
PCM240I	240	200 - 240	200 - 240	1.0	50/60	IEC320 / (4) IEC320	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	12 (5.4)
PCM360A	360	120	120	3.0	60	5-15P HG / (4) 5-15R HG	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	12 (5.4)
PCM360I	360	200 - 240	200 - 240	1.5	50/60	IEC320 / (4) IEC320	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	14.5 (6.6)
PCM360J	360	100 - 120	100 - 120	3.0	50/60	IEC320 / (4) IEC320	4.79 x 4.11 x 8.97 (12.2 x 10.4 x 22.8)	14.5 (6.6)
PCM550A	550	120	120	4.6	60	5-15P HG / (4) 5-15R HG	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	22 (10.0)
PCM550I	550	200 - 240	200 - 240	2.3	50/60	IEC320 / (4) IEC320	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	23.5 (10.7)
PCM750A	750	120	120	6.2	60	5-15P HG / (4) 5-15R HG	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	31 (14.1)

* Six foot line cord is standard on 60 Hz models. Additional configurations available (call factory).

** Also available with (6) 5-15R HG.

*** Input voltage range is selectable.

PCm Series Medical-Grade Power Conditioner

Model Number*	Load (VA)	Nominal Input (Volts)	Output (Volts)	Max. Load (Amps)	Frequency (Hz)	Connectors: Input / Output*	Dimensions (HxWxD) in (cm)	Ship Weight lbs (kg)
PCM750I	750	200 - 240	200 - 240	3.1	50/60	IEC320 / (4) IEC320	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	26.5 (12.0)
PCM750J	750	100-120	100-120	6.2	50/60	IEC320 / (4) IEC320	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	32 (14.5)
PCM1000A	1000	120	120	8.4	60	5-15P HG / (4) 5-15R HG**	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	31 (14.1)
PCM1000I	1000	200-240	200 - 240	4.2	50/60	IEC320 / (4) IEC320	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	31.5 (14.3)
PCM1000J	1000	100-120	100-120	8.4	50/60	IEC320 / (4) IEC320	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	32 (14.5)
PCM1000Y	1000	200-240	100-120	8.4	50/60	IEC320 / (4) IEC320	5.4 x 6.4 x 10.9 (13.7 x 16.4 x 27.7)	32 (14.5)
PCm1500A	1500	120	120	12.5	60	5-15P HG / (6) 5-15R HG	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm1500I	1500	200 - 240	200 - 240	6.3	50/60	IEC cord set / (6) IEC320	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm1500J	1500	***100 - 120 or 200 - 240	100 - 120	12.5	50/60	hardwired / hardwired	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm1500W	1500	***100 - 120 or 200 - 240	200 - 240	6.3	50/60	hardwired / hardwired	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm2000A	1920	120	120	16	60	5-20P HG / (4) 5-15R HG & (2) 5-20R HG	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)

* Six foot line cord is standard on 60 Hz models. Additional configurations available (call factory).

** Also available with (6) 5-15R HG.

*** Input voltage range is selectable.

CHLORIDE

PCm Series Medical-Grade Power Conditioner

Model Number*	Load (VA)	Nominal Input (Volts)	Output (Volts)	Max. Load (Amps)	Frequency (Hz)	Connectors: Input / Output*	Dimensions (HxWxD) in (cm)	Ship Weight lbs (kg)
PCm2000I	2000	200 -240	200 - 240	8.4	50/60	IEC cord set / (6) IEC320	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm2000Y	2000	200-240	100-200	16.7	50/60	IEC cord set/ (6) IEC320 (10A) & (1) IEC320 (16A)	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm2000U	2000	***100 - 120 or 200 - 240	100 - 120	16.7	50/60	hardwired / hardwired	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm2000W	2000	***100 - 120 or 200 - 240	200 - 240	8.4	50/60	hardwired / hardwired	8 x 7 x 16.3 (20 x 18 x 41.4)	53 (23)
PCm3000A	3000	120	120	25	60	L5-30P HG / L5-30R & (2) 5-20R & (2) 5-15R	8 x 7 x 16.3 (20 x 18 x 41.4)	65 (29)
PCm3000I	3000	200 -240	200 - 240	12.5	50/60	IEC cord set / (1) IEC320(16A) & (6) IEC320(10A)	8 x 7 x 16.3 (20 x 18 x 41.4)	65 (29)
Pcm3000U	3000	***100 - 120 or 200 - 240	100 - 120	25	50/60	hardwired / hardwired	8 x 7 x 16.3 (20 x 18 x 41.4)	65 (29)
PCm3000W	3000	***100 - 120 or 200 - 240	200 - 240	12.5	50/60	hardwired / hardwired	8 x 7 x 16.3 (20 x 18 x 41.4)	65 (29)

* Six foot line cord is standard on 60 Hz models. Additional configurations available (call factory).

** Also available with (6) 5-15R HG.

*** Input voltage range is selectable.